




| | |
|---|--|
|  | Rapid Competitive Immunochemical Test for detection of Methamphetamine in urine |
| Application | For the detection of presence of methamphetamine in human urine specimens qualitatively with concentrations as low as 1000 ng/ml |
| Principle | Self performing, two site sandwich immunoassay, lateral flow device |
| Reagents | <p>A. Each INSIGHT MET kit contains individual pouches each containing a</p> <ol style="list-style-type: none"> 1. Device: Membrane test assembly impregnated with colloidal gold conjugated to anti methamphetamine antibody monoclonal and mouse IgG, methamphetamine conjugated to BSA and goat anti mouse IgG at the respective regions. 2. Sample dropper. 3. Desiccant pouch. <p>B. Package insert.</p> |

| Presentation | Pack | Cat.No. |
|--|----------|----------|
| Insight™ MET kit contains individually pouched device | 10 Tests | 10804010 |

