



	Rapid Competitive Immunochromatographic Test for detection of methadone in urine
Application	For the detection of presence of methadone in human urine specimens qualitatively with concentrations as low as 300 ng/ml
Principle	Self performing, two site sandwich immunoassay, lateral flow device
Reagents	A. Each INSIGHT MTD kit contains individual pouches each containing a <ol style="list-style-type: none"> 1. Device: Membrane test assembly impregnated with colloidal gold conjugated to anti methadone antibody monoclonal and mouse IgG, methadone conjugated to BSA and goat anti mouse IgG at the respective regions. 2. Sample dropper. 3. Desiccant pouch. B. Package insert.

Presentation	Pack	Cat.No.
Insight™ MTD kit contains individually pouched device	10 Tests	10809010

