



	<b>Rapid Competitive Immunochemical Test for detection of Opiates in urine</b>
<b>Application</b>	For the detection of presence of morphine in human urine specimens qualitatively with concentrations as low as 300 ng/ml
<b>Principle</b>	Self performing, two site sandwich immunoassay, lateral flow device
<b>Reagents</b>	A. Each INSIGHT OPI kit contains individual pouches each containing a <ol style="list-style-type: none"> <li>1. Device: Membrane test assembly impregnated with colloidal gold conjugated to anti morphine antibody monoclonal and mouse IgG, morphine conjugated to BSA and goat anti mouse IgG at the respective regions.</li> <li>2. Sample dropper.</li> <li>3. Desiccant pouch.</li> </ol> B. Package insert.
<b>Test duration</b>	15 mins.
<b>Storage / Stability</b>	2-8°C, 18 months

Presentation	Pack	Cat.No.
<b>Insight™ OPI kit contains individually pouched device</b>	10 Tests	10802010

