



	Rapid Competitive Immunochemical Test for detection of Opiates in urine
Application	For the detection of presence of morphine in human urine specimens qualitatively with concentrations as low as 300 ng/ml
Principle	Self performing, two site sandwich immunoassay, lateral flow device
Reagents	<p>A. Each INSIGHT OPI kit contains individual pouches each containing a</p> <ol style="list-style-type: none"> 1. Device: Membrane test assembly impregnated with colloidal gold conjugated to anti morphine antibody monoclonal and mouse IgG, morphine conjugated to BSA and goat anti mouse IgG at the respective regions. 2. Sample dropper. 3. Desiccant pouch. <p>B. Package insert.</p>
Test duration	15 mins.
Storage / Stability	2-8°C, 18 months

Presentation	Pack	Cat.No.
Insight™ OPI kit contains individually pouched device	10 Tests	10801010