



	<b>Rapid Competitive Immuno-chromatographic Test for detection of Phencyclidine in urine</b>
<b>Application</b>	Detects the presence of Phencyclidine in human urine at concentrations as low as 25 ng/ ml.
<b>Principle</b>	Self performing, semi-quantitative two-site sandwich immunoassay, lateral flow device
<b>Reagents</b>	<p>A. Each INSIGHT PCP kit contains individual pouches each containing</p> <ol style="list-style-type: none"> <li>1. Device: Membrane test assembly impregnated with colloidal gold conjugated to anti-phencyclidine antibody monoclonal and rabbit IgG, phencyclidine conjugated to BSA and anti rabbit globulin at the respective regions.</li> <li>2. Sample dropper.</li> <li>3. Desiccant pouch.</li> </ol> <p>B. Package insert.</p>
<b>Test duration</b>	Five minutes assay
<b>Storage / Stability</b>	4-30°C, 12 months

Presentation	Pack	Cat.No.
<b>Insight™ PCP kit contains individually pouched device</b>	10 Tests	10811010

