

	<b>Rapid Competitive Immuno chromatographic Test for detection of Propoxyphene in urine</b>
<b>Application</b>	Detects the presence of <b>Propoxyphene</b> in human urine at concentrations as low as 300 ng/ ml.
<b>Principle</b>	Self performing, semi-quantitative two-site sandwich immunoassay, lateral flow device
<b>Reagents</b>	A. Each INSIGHT PPX kit contains individual pouches each containing <ol style="list-style-type: none"> <li>1. Device: Membrane test assembly impregnated with colloidal gold conjugated to anti-Propoxyphene antibody monoclonal and rabbit IgG, Propoxyphene conjugated to BSA and anti rabbit globulin at the respective regions.</li> <li>2. Sample dropper.</li> <li>3. Desiccant pouch.</li> </ol> B. Package insert.
<b>Test duration</b>	Five minutes assay
<b>Storage / Stability</b>	4-30°C, 12 months

Presentation	Pack	Cat.No.
<b>Insight™ PPX kit contains individually pouched device</b>	10 Tests	10812010

