



	Rapid Immunochromatographic assay for the detection of Rheumatoid factors in human serum
Application	For the detection of rheumatoid factors in human serum
Sensitivity	10 IU/ml
Reagents	<p>A. Each INSIGHT RF kit contains individual pouches each containing:</p> <ol style="list-style-type: none"> 1. Device: Membrane test assembly impregnated with colloidal gold conjugated to agglutinating sera for human IgM and rabbit globulin, Fc fraction of human IgG and agglutinating sera for rabbit globulin at the respective regions. 2. Desiccant pouch. 3. Sample applicator. 4. sample running buffer <p>B. Package insert.</p>
Test duration	Five minutes assay
Storage / Stability	4-30°C, 18 months

Presentation	Pack	Cat.No.
Insight™ RF kit contains individually pouched device	10 Tests	108180010
	25 Tests	108180025

