



	Rapid Competitive Immunochromatographic Test for detection of tetrahydrocannabinol in urine
Application	For the detection of presence of tetrahydrocannabinol in human urine specimens qualitatively with concentrations as low as 50 ng/ml
Principle	Self performing, two site sandwich immunoassay, lateral flow device
Reagents	<p>A. Each INSIGHT THC kit contains individual pouches each containing a</p> <ol style="list-style-type: none"> 1. Device: Membrane test assembly impregnated with colloidal gold conjugated to anti tetrahydrocannabinol antibody monoclonal and mouse IgG, tetrahydrocannabinol conjugated to BSA and goat anti mouse IgG at the respective regions. 2. Sample dropper. 3. Desiccant pouch. <p>B. Package insert.</p>

Presentation	Pack	Cat.No.
Insight™ THC kit contains individually pouched device	10 Tests	10807010

