



	Rapid Test for <i>P. falciparum</i> malaria
Application	For the detection of <i>P. falciparum</i> specific HRP 2 in whole blood
Principle	Self performing, two site sandwich immunoassay, lateral flow device
Sensitivity	> 99% correlation with microscopy
Specificity	100% specific to <i>P. falciparum</i> HRP 2
Reagents	<ul style="list-style-type: none"> • Membrane assembly predisposed with <ul style="list-style-type: none"> - Agglutinating Sera for HRP-2 - colloidal gold conjugat, rabbit IgG colloidal gold conjugate - Agglutinating Sera for HRP-2 at 'Test' region - Agglutinating Sera for rabbit globulin at 'Control' region • Clearing Buffer
Test duration	Twenty minutes assay
Storage / Stability	4-45°C, 2 years

Presentation	Pack	Cat.No.
PARACHECK Pf [®] kit contains individually pouched device with clearing buffer and 5 µl sample loop	1 Test	302030001
	5 Tests	302030005
	10 Tests	302030010
	25 Tests	302030025
	100 Tests	302030100