



Retrocheck HIV WB	Rapid immunochromatographic test for detection of HIV1/2 antibodies in serum, plasma and whole blood
Application	For simultaneous detection of antibodies specific to HIV 1/2 in serum, plasma and whole blood
Principle	<i>In vitro</i> , rapid, self performing, qualitative, two site, sandwich lateral flow immunochromatographic assay
Sensitivity	100% sensitive for HIV 1/2 antibodies
Specificity	99.75% specificity
Reagents	Membrane assembly comprises of <ul style="list-style-type: none"> • HIV 1/2 specific recombinant antigen- colloidal goldconjugate co-dispensed with rabbit IgG colloidal gold conjugate • HIV 1/2 specific recombinant antigen at region 'T' • Agglutinating sera for rabbit globulin at region 'C'.
Test duration	Fifteen minutes
Storage / Stability	4-30°C, 2 years

Presentation	Pack	Cat.No.
RETROCHECK HIV kit contains individually pouched devices and sample running buffer	1 TEST	402010001
	10 TESTS	402010010
	25 TESTS	402010025
	50 TESTS	402010050
	100 TESTS	402010100