



<p>Retrocheck HIV WB</p>	<p>Rapid immunochromatographic test for detection of HIV1/2 antibodies in serum, plasma and whole blood</p>
<p>Application</p>	<p>For simultaneous detection of antibodies specific to HIV 1/2 in serum, plasma and whole blood</p>
<p>Principle</p>	<p><i>In vitro</i>, rapid, self performing, qualitative, two site, sandwich lateral flow immunochromatographic assay</p>
<p>Sensitivity</p>	<p>100% sensitive for HIV 1/2 antibodies</p>
<p>Specificity</p>	<p>99.75% specificity</p>
<p>Reagents</p>	<p>Membrane assembly comprises of</p> <ul style="list-style-type: none"> • HIV 1/2 specific recombinant antigen- colloidal goldconjugate co-dispensed with rabbit IgG colloidal gold conjugate • HIV 1/2 specific recombinant antigen at region 'T' • Agglutinating sera for rabbit globulin at region 'C'.
<p>Test duration</p>	<p>Fifteen minutes</p>
<p>Storage / Stability</p>	<p>4-30°C, 2 years</p>

Presentation	Pack	Cat.No.
<p>RETROCHECK HIV kit contains individually pouched devices and sample running buffer</p>	<p>10 TESTS</p>	<p>40501010</p>
	<p>50 TESTS</p>	<p>402010050</p>
	<p>100 TESTS</p>	<p>402010100</p>



