



	Rapid Immunoconcentration Test for HIV1 and HIV 2 antibodies
Application	For detection of antibodies specific to HIV 1 and HIV 2 in human serum, plasma
Principle	Membrane based flow through immunoassay for the detection of HIV 1 and HIV 2 antibodies
Sensitivity	100% sensitive for HIV 1 and HIV 2 antibodies
Specificity	100% specific for HIV 1 and HIV 2 antibodies
Reagents	<ul style="list-style-type: none"> • Ready to use individually pouched, flow through test devices striped with • HIV 1 specific purified synthetic peptides at test region '1' • HIV 2 specific purified synthetic peptides at test region '2' • A blue dyed protein A based control band at region 'C'
Test duration	5 minutes
Storage / Stability	2-8°C, 12 months

Presentation	Pack	Cat.No.
RETROQUIC HIV kit contains individually pouched devices, wash buffer and protein A conjugated gold sol solution.	10 TESTS	40502010
	50 TESTS	402020050
	75 TESTS	40502075
	100 TESTS	402020100

