



One Step Test for HBsAg

Application	For the detection of Hepatitis B surface antigen, a marker for Hepatitis B infections in serum / plasma specimen.
Principle	Rapid, qualitative, two site sandwich immunoassay, lateral flow device
Detection limit	0.3 ng/ml of HBsAg
Sensitivity	100%
Specificity	100%
Reagents	Membrane assembly predisposed with: <ul style="list-style-type: none"> • Agglutinating sera for HBsAg colloidal gold conjugate, rabbit IgG colloidal gold conjugate • Agglutinating sera for HBsAg at 'Test' region • Agglutinating sera for rabbit globulin at 'Control' region
Test duration	20 mins.
Storage / Stability	4-30°C, 2 years

Presentation	Pack	Cat.No.
SENSAs Hep-B kit contains individually pouched device with sample dropper	25 TESTS	402150025