



SEROCHECK-Tp™	Rapid test for Syphilis
Application	Modified TPHA, for the detection of <i>Treponema</i> specific IgG and IgM antibodies in human serum or plasma
Principle	Self-performing, double antigen sandwich immunoassay, lateral flow dipstick
Sensitivity	100%
Specificity	98%
Reagents	<p>A. Dipstick container with individual pouches and desiccant. Each individual pouch contain Membrane assembly pre-dispensed with:</p> <ul style="list-style-type: none"> • Recombinant <i>Treponema pallidum</i> antigen colloidal gold conjugate. • Recombinant <i>Treponema pallidum</i> antigen at 'Test' region. • Agglutinating sera for rabbit globulin coated at 'Control' region. <p>B. Dipstick Holder.</p> <p>C. Package insert</p>
Test duration	Fifteen minutes assay
Storage / Stability	4-30°C

Presentation	Pack	Cat.No.
SEROCHECK Tp kit contains individually pouched dipsticks with desiccant.	25 Tests	501090025