



	<p>A rapid, qualitative, Immunochromatographic assay for the simultaneous detection of Hb S and Hb A in human whole blood sample.</p>
<p>Application</p>	<p>For simultaneous detection of Hb S and Hb A in human whole blood.</p>
<p>Sensitivity</p>	<p>97%</p>
<p>Specificity</p>	<p>Specific to Hb S and Hb A</p>
<p>Reagents</p>	<ul style="list-style-type: none"> • Sicklecheck kit contains: individual pouches each containing: <ol style="list-style-type: none"> 1. Device: Membrane test assembly impregnated with colloidal gold conjugated to Anti-Hb S, Anti-Hb A antibody, Hb S & Hb A and goat anti mouse IgG at the respective regions 2. Desiccant pouch. 3. Disposable plastic sample applicator. 4. Assay buffer • Package insert
<p>Test duration</p>	<p>Fifteen minutes assay</p>
<p>Storage / Stability</p>	<p>4-30°C, 24 months</p>

Presentation	Pack	Cat.No.
Sicklecheck contains individually pouched device	10 T	508010010