



<h1>Syphichек®</h1>	<h2>Rapid Test for Syphilis</h2>
Application	Modified TPHA, for the detection of <i>Treponema</i> specific IgG and IgM antibodies in serum or plasma
Principle	Self performing Syphilis, double antigen sandwich immunoassay, lateral flow dipstick
Sensitivity	100% correlation with TPHA
Specificity	Recombinant <i>Treponema</i> antigens, specific for IgG and IgM class of Anti - treponemal antibodies
Reagents	Membrane assembly pre-dispensed with <ul style="list-style-type: none"> • Recombinant <i>T. pallidum</i> antigen - colloidal gold conjugate • Recombinant <i>T. pallidum</i> antigen at 'Test' region • Agglutinating sera for rabbit globulin at 'Control' region • Clearing Buffer
Test duration	Fifteen minutes assay
Storage / Stability	4-30°C, 2 years

Presentation	Pack	Cat.No.
Syphichек® kit contains individually pouched dipstick with sample dropper and clearing buffer	25 Tests	401020025