



Syphicheck-WB[®]	Rapid Test for Syphilis
Application	Modified TPHA, for the detection of <i>Treponema</i> specific IgG and IgM antibodies in whole blood, serum or plasma
Principle	Self performing, double antigen sandwich immunoassay, lateral flow device
Sensitivity	100% correlation with TPHA
Specificity	Recombinant <i>Treponema</i> antigens, specific for IgG and IgM class of Anti - treponemal antibodies
Reagents	Membrane assembly comprises of <ul style="list-style-type: none"> • Recombinant <i>T. pallidum</i> antigen - colloidal gold conjugate • Recombinant <i>T. pallidum</i> antigen at 'Test' region • Agglutinating sera for rabbit globulin at 'Control' region • Clearing Buffer
Test duration	Fifteen minutes assay
Storage / Stability	4-30°C, 2 years

Presentation	Pack	Cat.No.
SYPHICHECK-WB[®] kit contains individually pouched device with sample dropper and clearing buffer	25 Tests	401030025
	100 Tests	401030100

