



	<b>Rapid One Step Test (dipstick) for detection of HCG</b>
<b>Application</b>	For the detection of human chorionic gonadotrophin (HCG) hormone in urine
<b>Principle</b>	Self performing, two site sandwich immunoassay, lateral flow dipstick
<b>Sensitivity</b>	~ 25 mIU/ml
<b>Reagents</b>	Each kit contains: A. Individual pouches, each containing a <ol style="list-style-type: none"> <li>1. Membrane assembly impregnated with Anti HCG antiserum- colloidal gold conjugate, rabbit IgG colloidal gold conjugate, anti-HCG antiserum and anti rabbit antiserum coated at the respective regions.</li> <li>2. Disposable Plastic Dropper.</li> <li>3. Desiccant pouch.</li> </ol> B. Package insert.
<b>Test duration</b>	5 mins.
<b>Storage / Stability</b>	4-30°C, 2 years

Presentation	Pack	Cat.No.
<b>Tina™ HCG</b> kit contains individually pouched dipstick	30 Tests	103300030

