



<b>Toxocheck™</b>	<b>Rapid immunochromatographic test for IgM and IgG antibodies to <i>T. gondii</i></b>
<b>Application</b>	Rapid test for detection of IgM & IgG antibodies to <i>Toxoplasma gondii</i> in human serum or plasma.
<b>Principle</b>	Rapid, qualitative, two site sandwich immunoassay, lateral flow device
<b>Sensitivity</b>	100%
<b>Specificity</b>	98.58%
<b>Reagents</b>	Membrane assembly for <i>Toxoplasma</i> IgG/IgM antigen is predisposed with: <ul style="list-style-type: none"> <li>• Recombinant <i>T. gondii</i> antigen–colloidal gold conjugate and streptavidin gold conjugate</li> <li>• Anti-human IgG and Anti-human IgM antibodies are immobilized as two individual test bands at region 'G' and region 'M' respectively</li> <li>• Biotinylated BSA predisposed at the control window 'C'</li> </ul>
<b>Test duration</b>	15 mins
<b>Storage / Stability</b>	4-30°C, 24 months

<b>Presentation</b>	<b>Pack</b>	<b>Cat.No.</b>
<b>Toxocheck™ I</b> kit contains individually pouched device and sample running buffer	10 Tests	503170010
	25 Tests	503170025

