



	<h2>Rapid One Step Test for HBsAg</h2>
Application	For the detection of HBsAg in serum / plasma
Principle	Self performing, two site sandwich immunoassay, lateral flow dipstick
Sensitivity	~ 0.5 ng/ml of HBsAg
Specificity	100% specific to HBsAg (ad/ay subtypes)
Reagents	Membrane assembly predisposed with: <ul style="list-style-type: none"> - Agglutinating sera for HBsAg - colloidal gold conjugate, rabbit IgG colloidal gold conjugate - Agglutinating sera for HBsAg at 'Test' region - Agglutinating sera for rabbit globulin at 'Control' region
Test duration	15 mins.
Storage / Stability	4-30°C, 2 years

Presentation	Pack	Cat.No.
VIRUCHECK [®] kit contains individually pouched dipstick	1 Test	302020001
	10 Tests	302020010
	25 Tests	302020025
	50 Tests	302020050
	100 Tests	302020100