

**Brucel<sup>®</sup> - A**Test for the detection of *Brucella abortus* antibodies in slide and tube test

Intended Use: \_\_\_\_\_

Brucel<sup>®</sup> - A, a test for the detection of *Brucella abortus* antibodies in slide and tube test.Brucel<sup>®</sup> - A Kit components:

Reagent	Ready to use standardized, attenuated, stained, smooth specific antigen suspensions of <i>Brucella abortus</i>
Other Accessories	Package Insert

Performance Characteristics:

Specificity: Specific reactivity towards *Brucella abortus*

FEATURES	BENEFITS
Thermostabilizer - based reagents	New revolutionary thermostabilizer - based reagents makes them robust ensuring better kit performance under varying ambient conditions
Standardized, specific and smooth antigens	Ensures maximum sensitivity & specificity
Stringent GMP - based QC	Ensures lot-to-lot consistency in performance
Suitable for both slide and tube test	Provides versatility to the reagents
Long shelf life	Prevents reagent wastage : reducing costs

Storage/Stability	Temperature	Duration
Unopened kit	2-8°C	18 Months
Opened kit	2-8°C	18 Months

Available Pack Sizes
5 ml