



Brucel[®] - RB
Slide screening test for *Brucella* antibodies

Intended Use: _____

Brucel[®] - RB, a test for the detection of antibodies to *Brucella* antigens.

Brucel[®] - RB Kit components:

Reagent	Smooth, killed buffered suspensions of <i>Brucella abortus</i> strain 99, colored with rose bengal
Other Accessories	Package Insert

Performance Characteristics:

Specificity: Specific reactivity towards antibodies to *Brucella* antigens

FEATURES	BENEFITS
Uses Internationally recommended <i>B. abortus</i> strain - 99	High sensitivity and specificity
Rose Bengal stain standardized with 2 nd International preparation	Reliable performance

Storage/Stability	Temperature	Duration
Unopened kit	2-8°C	18 Months
Opened kit	2-8°C	18 Months

Available Pack Sizes
5 ml