

**Foretel: Latex agglutination Inhibition Slide Test for hCG****Intended Use:** _____

Foretel, a latex agglutination inhibition slide test for hCG which is an indirect test for the detection of pregnancy in an urine sample.

Foretel kit components:

Latex Reagent	Blend of small and large polystyrene ULP (uniform latex particles) coated with purified hCG
Antibody Reagent	beta hCG specific IgG class of murine monoclonal antibodies
Slides	Glass slides with three reaction circles
Other Accessories	Droppers, mixing sticks, rubber teats and Package Insert

Performance Characteristics:

Sensitivity: ~ 0.3 IU/ml of hCG	Specificity: 100% specific to beta hCG
--	---

FEATURES	BENEFITS
Sensitivity 0.3 IU/ml	Quick detection
Latex reagent consists of a blend of small and large polystyrene uniform latex particles coated with purified hCG	Unique combination of small and large Uniform latex particles improves the avidity and reliability of the test
Antibody reagent consists of beta hCG specific murine monoclonal (IgG)	High specificity, eliminates cross reactivity with LH, FSH and TSH hormones
Distinct clear cut agglutination	Better interpretation of results

Storage / Stability	Temperature	Duration
Unopened kit	2–8°C	24 Months
Opened kit	2–8°C	24 Months

Available Pack Sizes		
25 Tests	35 Tests	70 Tests
100 Tests	140 Tests	