



## Paratrol™ Controls for validation of Malaria Rapid Tests

**Intended Use:** \_\_\_\_\_

Controls for on site validation / quality assurance of immunological assays for malaria antigen detection.

**Paratrol™ Kit components:**

<b>Reagent</b>	Set of lyophilized preparation of human whole blood sources parasitized with <i>P.falciparum</i> , <i>P. vivax</i> parasites and malaria parasite negative whole blood
<b>Other Accessories</b>	Package Insert

FEATURES	BENEFITS
<b>Precisely prepared reagent-with known parasitic load</b>	Appropriately reconstituted controls, ensures quality and standardized results of malaria antigen detection immuno assay
<b>Freeze dried reagent</b>	Aids in minimizing microbial contamination
<b>Convenient pack size</b>	Ideal for all users

Storage / Stability	Temperature	Duration
<b>Unopened kit</b>	2-8°C	6 Months
<b>Reconstituted vial</b>	2-8°C	7 Days

Available Pack Sizes
3 x 2 x 0.2 ml