



Rhelax[®] ASO

Latex agglutination slide test for detection of Anti-Streptolysin O

Intended Use: _____

Rhelax[®] ASO, a latex agglutination slide test for the detection and semiquantitation of Anti - Streptolysin O in serum.

Rhelax[®] ASO Kit components:

Latex Reagent	Uniform suspension of polystyrene latex particles coated with stabilized Streptolysin O
Positive Control	Reactive with the reagent
Negative Control	Non Reactive with the reagent
Slide	Slide with six reaction circles
Other Accessories	Sample droppers, mixing sticks, rubber teat and Package Insert

Performance Characteristics:

Sensitivity: ≤ 200 IU/ml	Specificity: Purified Streptolysin O specific for Anti-Streptolysin O antibodies
---------------------------------	---

FEATURES	BENEFITS
Sensitivity accurately calibrated at 200 IU/ml	Excludes reactivity with normal samples Reacts with true positives
Prozone limit > 4000 IU/ml	Wide working range No risk of false negatives during single screens Confident reporting of results
Stabilized SLO used on latex	Improved in - use reagent stability and sensitivity
Enhanced readability	Smooth negatives Crisp positives

Storage / Stability	Temperature	Duration
Unopened kit	2-8°C	18 Months
Opened kit	2-8°C	18 Months

Available Pack Sizes		
5 ml	25 Tests	35 Tests
50 Tests	70 Tests	100 Tests