



Rhelax[®] CRP

Latex agglutination slide test for detection of C - Reactive Protein

Intended Use: _____

Rhelax[®] CRP, a latex agglutination slide test for the detection and semiquantitation of C - Reactive Protein in serum.

Rhelax[®] CRP Kit components:

Latex Reagent	Uniform suspension of polystyrene latex particles coated with Anti-CRP antibodies
Positive Control	Reactive with the reagent
Negative Control	Non Reactive with the reagent
Slide	Slide with six reaction circles
Other Accessories	Sample droppers, mixing sticks, rubber teat and Package Insert

Performance Characteristics:

Sensitivity: ~ 0.6 mg/dl	Specificity: Anti-CRP antibodies of IgG class specific to CRP
---------------------------------	--

FEATURES	BENEFITS
Accurate calibration at 0.6 mg/dl	Excludes reactivity with normal samples No false positives Correct classification of titres
High capacity and avid antibodies	No prozone till 100 mg/dl Wide working range Reliable detection
Enhanced readability	Smooth negatives Crisp positives

Storage / Stability	Temperature	Duration
Unopened kit	2-8°C	18 Months
Opened kit	2-8°C	18 Months

Available Pack Sizes		
5 ml	25 Tests	35 Tests
50 Tests	70 Tests	100 Tests