



## Rhelax® SLE

### Latex agglutination slide test for Anti-DNP in Systemic Lupus Erythematosus

**Intended Use:** \_\_\_\_\_

Rhelax® SLE, a latex agglutination slide test for the detection of serum Anti-DNP (Anti-Deoxyribonucleoprotein) in Systemic Lupus Erythematosus (SLE).

**Rhelax® SLE Kit components:**

<b>Latex Reagent</b>	Uniform suspension of polystyrene latex particles coated with Deoxyribonucleoprotein
<b>Positive Control</b>	Reactive with the reagent
<b>Negative Control</b>	Non Reactive with the reagent
<b>Slide</b>	Slide with six reaction circles
<b>Other Accessories</b>	Sample droppers, mixing sticks, rubber teat and Package Insert

**Performance Characteristics:**

<b>Sensitivity:</b> 100 %	<b>Specificity:</b> Deoxyribonucleoprotein specific for Anti-DNP (DNA + Histone complex) antibodies
---------------------------	---

FEATURES	BENEFITS
<b>Specificity: Anti DNP</b>	Reactivity against Anti DNA and Anti histone antibodies (Anti DNP)
<b>Sensitivity</b>	Picks up Active SLE and drug induced SLE positivity
<b>Three minute test</b>	Rapid screening
<b>Complete kit</b>	Complete system for easy screening with positive and negative controls
<b>Pack size</b>	Convenient test pack with 18 months shelf life

Storage / Stability	Temperature	Duration
Unopened kit	2-8°C	18 Months
Opened kit	2-8°C	18 Months

Available Pack Sizes	
10 Tests	25 Tests