

Urinary Congo Red Spot Test for Pre-eclampsia (PE) Screening

Intended Use: _____

REDSPOT-PE™ is a point-of-care urine test, to detect urinary congophilia useful for screening of Pre-eclampsia.

REDSPOT-PE™ kit components:

Reagent	Ready to use Congo Red Reagent
Other Accessories	Test Cards, Sample-Reagent Mixing Vials, Sample Droppers, Test Droppers, Pack insert

Performance Characteristics:

Specificity: The test is based on the chemical interaction between Congo red (CR) and amyloid proteins (misfolded proteins) from urine sample.
--

Features	Benefits
Point of Care Urine test	Non-invasive, Easy to use, suitable for labs and OPD settings
Easy to Use with Visual Output	Minimal training required; test results are simple to interpret visually

Storage / Stability	Temperature	Duration
Unopened kit	25 – 30°C	12 months

Available Pack Sizes
50 Tests