

**Tulip's Febrile Antigen Set: Febrile Antigen Panel****Intended Use:** _____

Tulip's Febrile Antigen Set, a Febrile Antigen Panel for the detection of antibodies in certain febrile diseases such as Salmonellosis, Brucellosis and Rickettsial diseases in serum.

Tulip's Febrile Antigen Set kit components:

Reagent	Ready to use, colour coded, attenuated suspension of <i>S. typhi</i> 'O', <i>S. typhi</i> 'H', <i>S. paratyphi</i> 'AH', <i>S. paratyphi</i> 'BH', <i>Brucella abortus</i> , Proteus OX19 antigens for slide and tube tests
Positive Control	Polyspecific, reactive with <i>S. typhi</i> , <i>S. paratyphi</i> , <i>Brucella abortus</i> and Proteus OX19 antigens
Negative Control	Febrile antigen negative control
Other Accessories	Package Insert

Performance Characteristics:

Specificity: Detects antibodies in certain febrile diseases such as Salmonellosis, Brucellosis and Rickettsial diseases in serum

FEATURES	BENEFITS
All in one kit	Detects and differentiates the three febrile diseases using one kit
Convenient Pack sizes	No wastage

Storage / Stability	Temperature	Duration
Unopened kit	2–8°C	18 Months
Opened kit	2–8°C	18 Months

Available Pack Sizes
6 X 5 ml