



Typhocheck® Widal antigen for tube test

Intended Use: _____

Typhocheck® is Widal antigen for tube test which is used for the detection of antibodies to *S. typhi* 'O', *S. typhi* 'H', *S. paratyphi* 'AH', *S. paratyphi* 'BH', *S. paratyphi* 'CH', *S. paratyphi* 'AO', *S. paratyphi* 'BO', *S. paratyphi* 'CO' in serum.

Typhocheck® Kit components:

Reagent	Ready to use, color coded, attenuated suspension of <i>S. typhi</i> 'O', <i>S. typhi</i> 'H', <i>S. paratyphi</i> 'AH', <i>S. paratyphi</i> 'BH', <i>S. paratyphi</i> 'CH', <i>S. paratyphi</i> 'AO', <i>S. paratyphi</i> 'BO', <i>S. paratyphi</i> 'CO' antigens for tube tests
Other Accessories	Package Insert

Performance Characteristics:

Specificity: For the detection of antibodies to *S. typhi* and *S. paratyphi*

FEATURES	BENEFITS
Thermostabilizer - based reagents	New revolutionary thermostabilizer - based reagents makes them robust ensuring better kit performance under varying ambient conditions
Standardized, specific and smooth antigens	Ensures maximum sensitivity & specificity
Specialized staining of antigen as per International recommendation	Imparts high readability
Differential color coding of antigen vials / labels	Prevents performance errors
Multiple pack sizes	Offers choice of application
Polyspecific positive control provided with the kit	Can perform quality control & reagent validation

Storage/Stability	Temperature	Duration
Unopened kit	2-8°C	12 Months
Opened kit	2-8°C	12 Months

Available Pack Sizes
50 ml each
4 x 50 ml set