

**Turbodyne™ Control****Turbodyne™ Controls Kit components:**

Reagent	Lyophilised /Liquid stable controls prepared from Human serum/plasma/whole blood fortified with preservatives and stabilizers
Other Accessories	Package Insert

PRODUCT	ANALYTE CONTENT	Available Pack Size
Turbodyne™ RF Control	Rheumatoid Factor	0.5 ml
Turbodyne™ CRP Control	C-reactive protein	0.5 ml
Turbodyne™ ASO Control	Antistreptolysin O	0.5 ml
Turbodyne™ MA Control	Microalbumin	0.5 ml
Turbodyne™ D-Dimer Control	D-dimer	0.5 ml
Turbodyne™ Cystatin C Control	Cystatin C	0.2 ml
Turbodyne™ HbA1c Control	HbA1c	0.2 ml
Turbodyne™ IgE Control	Immunoglobulin-E	0.2 ml

FEATURES	BENEFITS
Convenient Pack size	Ideal for all users
Precisely Prepared Reagent	When reconstituted appropriately, ensures precise results

Storage / Stability	Temperature	Duration
Unopened kit	2-8°C	12 months
Lyophilised control after Reconstitution	2-8°C	1 week
Liquid stable Control	2-8°C	1 week