



Vital Widal[®] Vitaly stained Widal antigen for tube test

Intended Use: _____

Vital Widal[®] is vitally stained Widal antigen for tube test which is used for the detection of antibodies to *S. typhi* 'O', *S. typhi* 'H', *S. paratyphi* 'AH', *S. paratyphi* 'BH' in serum.

Vital Widal[®] Kit components:

Reagent	Ready to use, vitally stained smooth suspension of <i>S. typhi</i> 'O', <i>S. typhi</i> 'H', <i>S. paratyphi</i> 'AH', <i>S. paratyphi</i> 'BH' antigens for tube test
Positive Control	Polyspecific, reactive with <i>S. typhi</i> and <i>S. paratyphi</i> antigens
Other Accessories	Package Insert

Performance Characteristics:

Specificity: For the detection of antibodies to <i>S. typhi</i> and <i>S. paratyphi</i>

FEATURES	BENEFITS
Stain- while- they-grow Technology	Antigens are stained while they grow (Vital staining). The stain is taken up internally by bacterial cell
Vitaly stained antigens ensure no background colour	Clear-cut test results in tube test
Ideal as a confirmatory test for typhoid infection	Better than any other tube test

Storage / Stability	Temperature	Duration
Unopened kit	2-8°C	12 Months
Opened kit	2-8°C	12 Months

Available Pack Sizes	
4 x10 ml	2 x10 ml