



## BRUCEL-A

### Slide and tube test for *Brucella* antibodies

<b>Application</b>	For the detection of antibodies to <i>Brucella abortus</i>
<b>Specificity</b>	Specific reactivity towards <i>Brucella abortus</i> (Brucel-A)
<b>Reagent</b>	Ready to use, standardized, attenuated, stained, smooth specific antigen suspensions of <i>Brucella</i>
<b>Positive Control</b>	Reactive with the reagent
<b>Test duration</b>	<ul style="list-style-type: none"> <li>• Slide Test - one minute assay</li> <li>• Tube Test - 24 hours (overnight)</li> </ul>
<b>Storage / Stability</b>	2-8°C, 24 months

Presentation	Pack	Cat.No.
BRUCEL-A	5 ml	105610005

