



BRUCEL-M	Slide and tube test for Brucella antibodies
Application	For the detection of antibodies to <i>Brucella melitensis</i>
Specificity	Specific reactivity towards <i>Brucella melitensis</i> (Brucel-M)
Reagent	Ready to use, standardized, attenuated, stained, smooth specific antigen suspensions of Brucella
Positive Control	Reactive with the reagent
Test duration	<ul style="list-style-type: none"> • Slide Test - one minute assay • Tube Test - 24 hours (overnight)
Storage / Stability	2-8°C, 24 months

Presentation	Pack	Cat.No.
BRUCEL-M	5 ml	105620005

