



	<h2>RPR Card Test</h2>
Application	For the detection of Anti-Lipoidal antibodies in serum and plasma during Syphilis
Specificity	Lipid complex reactive with Anti-Lipoidal antibodies
Reagent	Microfine particulate carbon suspension coated with lipid complex
Positive Control	Reactive with the reagent
Negative Control	Non reactive with the reagent
Test duration	Eight minutes assay
Storage / Stability	2-8°C, 18 months

Presentation	Pack	Cat.No.
CARBOGEN® kit contains reagent and accessories *CARBOGEN® 5 ml reduced pack contains only reagent	5 ml*	10514005
	50 Tests	10514050
	100 Tests	10514100
	250 Tests	10514250
RPR Positive control	0.5 ml	105160005
	1.0 ml	105160001

