



## Latex agglutination Inhibition Slide Test for hCG

<b>Application</b>	Indirect test for the detection of pregnancy in urine
<b>Sensitivity</b>	~0.3 IU/ml of hCG
<b>Specificity</b>	100% specific to b hCG
<b>Latex reagent</b>	Blend of small and large polystyrene ULP (uniform latex particles) coated with purified hCG
<b>Antibody reagent</b>	b hCG specific IgG class of murine monoclonal antibodies
<b>Test duration</b>	Three minutes assay
<b>Storage / Stability</b>	2-8°C, 24 months

<b>Presentation</b>	<b>Pack</b>	<b>Cat.No.</b>
<b>FORETEL kit contains reagents and accessories</b>	25 Tests	10310025
	35 Tests	10310035
	50 Tests	10310050
	70 Tests	10310070
	100 Tests	10310100
	140 Tests	10310140

