



Haemagglutination Slide Test

Application	For the detection of heterophile antibodies in serum during Infectious Mononucleosis
Specificity	Specially treated horse erythrocytes with high specificity for heterophile antibodies associated with Infectious Mononucleosis
Reagent	Stabilized suspension of horse erythrocytes
Positive Control	Reactive with the reagent
Negative Control	Non reactive with the reagent
Test duration	Two minutes assay
Storage / Stability	2-8°C, 24 months

Presentation	Pack	Cat.No.
IMMUTEX kit contains reagents and accessories	20 Tests	10810020

