



	<b>Controls for validation of Malaria Rapid Tests</b>
<b>Application</b>	Controls for on site validation / quality assurance of immunological assays for malaria antigen detection.
<b>Reagents</b>	Set of lyophilized preparation of human whole blood sources parasitized with <i>P. falciparum</i> , <i>P. vivax</i> parasites and malaria parasite negative whole blood.
<b>Reagents</b>	<ul style="list-style-type: none"> <li>• Pf Positive Control : Pf HRP 2, Pan pLDH, Aldolase.</li> <li>• Pv Positive Control : Pv pLDH, Pan pLDH, Aldolase.</li> <li>• Malaria Negative Control : No malaria specific antigen</li> </ul>
<b>Storage / Stability</b>	2-8°C, 6 months

Presentation	Pack	Cat.No.
PARATROL	3x2x0.2ml	110456002

