



	<b>Polyspecific Positive control</b>
<b>Application</b>	For validation of the performance of Progen OX19, OXK and OX2 antigen suspensions.
<b>Specificity</b>	Specific reactivity towards Progen OX19, OXK and OX2
<b>Reagent</b>	Ready to use standardised Goat antiserum with polyspecific antibodies
<b>Storage / Stability</b>	2-8°C, 18 months

Presentation	Pack	Cat.No.
<b>PROGEN polyspecific positive control</b>	0.5 ml	105840005
	1 ml	105840001

