



PUO SCREEN

Screening test for Typhoid, Brucellosis and Scrub typhus fever

Application	For the detection of antibodies to <i>S. typhi</i> 'O', <i>Brucella abortus</i> and <i>Proteus OXK</i> antigen in serum
Reagent	Ready to use, colour coded, attenuated suspension of <i>S. typhi</i> 'O', <i>Brucella abortus</i> and <i>Proteus OXK</i> antigen for slide test
Positive Control	Polyspecific, reactive with <i>S. typhi</i> 'O', <i>Brucella abortus</i> and <i>Proteus OXK</i> antigens
Test duration	Slide Test - one minute assay
Storage / Stability	2-8°C, 18 months

Presentation	Pack	Cat.No.
PUO Screen Set contains reagents and accessories	3 x 2 ml	105910032