



	Proteus antigen suspensions for Weil-Felix Test
Application	For the diagnosis of Rickettsial infection and differential diagnosis in patients with febrile fever.
Specificity	Specific reactivity towards Rickettsial antibodies.
Reagent	Ready to use, standardized, killed, stained, smooth specific antigen suspensions of Proteus OXK
Test duration	<ul style="list-style-type: none"> • Slide Test - one minute assay • Tube Test - 18 hours (overnight)
Storage / Stability	2-8°C, 18 months

Presentation	Pack	Cat.No.
PROGEN OX 2	5 ml	105830005

