



	<h2>Latex Agglutination Slide Test</h2>
Application	For the detection and semiquantitation of Anti-Streptolysin O in serum
Sensitivity	≥ 200 IU/ml
Specificity	Purified Streptolysin O specific for Anti-Streptolysin O antibodies
Latex reagent	Uniform suspension of polystyrene latex particles coated with stabilized Streptolysin O
Positive Control	Reactive with the reagent
Negative Control	Non reactive with the reagent
Test duration	Two minutes assay
Storage / Stability	2-8°C, 18 months

Presentation	Pack	Cat.No.
RHELAX-ASO kit contains reagents and accessories	5 ml	10430005
	25 Tests	10430025
	35 Tests	104300035
	50 Tests	10430050
	70 Tests	104300070
	100 Tests	10430100