



## Latex Agglutination Slide Test

<b>Application</b>	For the detection and semiquantitation of IgM rheumatoid factors in serum
<b>Sensitivity</b>	~ 10 IU/ml
<b>Specificity</b>	Suitably modified Fc fraction of human IgG specific for IgM rheumatoid factors
<b>Latex reagent</b>	Uniform suspension of polystyrene latex particles coated with suitably modified Fc fraction of IgG
<b>Positive Control</b>	Reactive with the reagent
<b>Negative Control</b>	Non reactive with the reagent
<b>Test duration</b>	Two minutes assay
<b>Storage / Stability</b>	2-8°C, 18 months

Presentation	Pack	Cat.No.
<b>RHELAX-RF kit contains reagents and accessories</b>	5 ml	10410005
	25 Tests	10410025
	35 Tests	104100035
	50 Tests	10410050
	70 Tests	104100070
	100 Tests	10410100