



## Latex Agglutination Slide Test

<b>Application</b>	For the detection of serum Anti-DNP (Anti-Deoxyribonucleoprotein) in Systemic Lupus Erythematosus (SLE)
<b>Specificity</b>	Deoxyribonucleoprotein specific for Anti-DNP (DNA + Histone complex) antibodies
<b>Latex reagent</b>	Uniform suspension of polystyrene latex particles coated with Deoxyribonucleoprotein
<b>Positive Control</b>	Reactive with the reagent
<b>Negative Control</b>	Non reactive with the reagent
<b>Test duration</b>	Two minutes assay
<b>Storage / Stability</b>	2-8°C, 18 months

Presentation	Pack	Cat.No.
<b>RHELAX-SLE® kit contains reagents and accessories</b>	10 Tests	10440010
	25 Tests	10440025
	50 Tests	10440050

