



REDSPOT-PE™	Point of care Urine test
Application	Urinary Congo red Test for Pre-eclampsia (PE) Screening.
Specificity	The test is based on the chemical interaction between Congo red (CR) and amyloid proteins (misfolded proteins) from urine sample with Cellulose from test card.
Reagent	<ul style="list-style-type: none"> • Ready to use Congo red reagent. • Test Cards (TC) with cellulose test paper.
Other accessories	Sample-Reagent Mixing Vials, Sample Droppers, Test Droppers, Pack insert
Test duration	10 minutes assay
Storage / Stability	25-30°C, 12 months

Presentation	Pack	Cat.No.
REDSPOT-PE™	50 Test	1108230050

