


**STAPHY LATEX TEST**
**Latex Agglutination Slide Test**

<b>Application</b>	For the Presumptive Identification of Staphylococcus Aureus
<b>Sensitivity</b>	95.8%
<b>Specificity</b>	100%
<b>Latex reagent</b>	Uniform suspension of polystyrene latex particles coated with human IgG and fibrinogen
<b>Positive Control</b>	Reactive with the latex test reagent and Non-reactive with the control reagent
<b>Control reagent</b>	Uncoated, plain, uniform suspension of polystyrene latex particles
<b>Test duration</b>	Two minutes assay
<b>Storage / Stability</b>	2-8°C, 24 months

<b>Presentation</b>	<b>Pack</b>	<b>Cat.No.</b>
<b>Staphy Latex kit contains reagents and accessories</b>	50 Tests	104500050