



	Latex Agglutination Slide Test
Application	For the screening and semiquantitation of antibodies to <i>T.gondii</i> in serum
Sensitivity	4 IU/ml
Specificity	<i>T. gondii</i> antigen specific to Anti-Toxoplasma antibodies of IgG and IgM class
Latex reagent	Uniform suspension of polystyrene latex particles coated with <i>T.gondii</i> antigen
Positive Control	Reactive with the reagent
Negative Control	Non reactive with the reagent
Test duration	Five minutes assay
Storage / Stability	2-8°C, 12 months

Presentation	Pack	Cat.No.
TOXOGEN® kit contains reagents and accessories	10 Tests	10570005
	25 Tests	10570025
	50 Tests	10570050

