



Modified VDRL Reagent

Application	For the detection of Anti-Lipoidal antibodies in serum, plasma and CSF during Syphilis
Specificity	Stabilized complex of cardiolipin, lecithin and cholesterol reactive with Anti-Lipoidal antibodies
Reagent	Ready to use stabilized emulsion of cardiolipin, lecithin and cholesterol
Positive Control	Reactive with the reagent
Negative Control	Non reactive with the reagent
Test duration	Four minutes assay
Storage / Stability	2-8°C, 18 months

Presentation	Pack	Cat.No.
TREPOLIPIN® kit contains reagent and controls	50 Tests	10511050
	100 Tests	10511100
	5 x 5 ml	10511505

