



## Modified VDRL Reagent

<b>Application</b>	For the detection of Anti-Lipoidal antibodies in serum, plasma and CSF during Syphilis
<b>Specificity</b>	Stabilized complex of cardiolipin, lecithin and cholesterol reactive with Anti-Lipoidal antibodies
<b>Reagent</b>	Ready to use stabilized emulsion of cardiolipin, lecithin and cholesterol
<b>Positive Control</b>	Reactive with the reagent
<b>Negative Control</b>	Non reactive with the reagent
<b>Test duration</b>	Four minutes assay
<b>Storage / Stability</b>	2-8°C, 12 months

Presentation	Pack	Cat.No.
<b>TREPOLIPIN® kit contains reagent and controls</b>	50 Tests	10511050
	100 Tests	10511100
	5 x 5 ml	10511505