



	<b>Febrile antigen panel</b>
<b>Application</b>	For the detection of antibodies in certain febrile diseases such as salmonellosis, brucellosis and rickettsial diseases in serum
<b>Reagent</b>	Ready to use, colour coded, attenuated suspension of <i>S. typhi</i> 'O', <i>S. typhi</i> 'O', <i>S. typhi</i> 'H', <i>S. paratyphi</i> 'AH', <i>S. paratyphi</i> 'BH', <i>Brucella abortus</i> , <i>Proteus OX19</i> antigens for slide and tube tests
<b>Positive Control</b>	Polyspecific, reactive with <i>S. typhi</i> , <i>S. paratyphi</i> , <i>Brucella abortus</i> and <i>Proteus OX19</i> antigens
<b>Test duration</b>	<ul style="list-style-type: none"> <li>• Slide Test - one minute assay</li> <li>• Tube Test - 18 to 20 hours (overnight)</li> </ul>
<b>Storage / Stability</b>	2-8°C, 18 months

Presentation	Pack	Cat.No.
<b>TULIP's FEBRILE ANTIGEN</b> Set contains reagents and accessories	6 x 5 ml	10591065

